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New Product JP PAT.No.5055524

This soil regulator is newly designed with totally new concept.

Main features

- Crushing and the sieving of the soil of φ 2mm or less can be done at the same time.
- All the processes can be done without putting out dust at all. Therefore, there is no worry of the health hazard by inhaling dust.
- There is no soil dust go out from equipment, so the dust collector is unnecessary.
- Since this equipment is compact design, it is possible to set up them in the most of laboratory.
- Because all the process are done in the airtight container, the cleaning time will be greatly reduced. It is used to taking long time to clean up soil dust more than the crushing and the sieving of the soil samples.
- Working hours can be shortened, and it correspond to a large amount of sample processing by improving the working efficiency of crushing and the riddle division process.

Specifications

| | |
|---------------------|--|
| Outside dimensions | W750×D670×H700mm (Exclude projection) |
| Weight | Approx. 160kg |
| Operation time | 30～120sec |
| Adjustable speed | 10～150rpm/sec |
| Rotational speed | 100～1000rpm |
| Angle of gradient | 0～30° |
| Power supply | AC100V 1A |
| Processing capacity | Air dried soil 150g / container × 4 = 600g |

Component

| | |
|--|-------|
| Dust shielded Automated mill and screen for soil | 1 |
| Exclusive container (Air tight container) | 8 set |
| Crushing rod (L50mm × φ18mm) | 8 |

Outer view

Inner view

< The condition of soil crushing >

Before processing

Air dried red grain soil (About10mm diameter) will be processed 90sec.

After processing

After the 90sec processing, all crushed fine soil particles shifted through 2mm screen. (Each soil particle diameter is Less than 2mm)

< Exclusive container >

Completely assembled exclusive container

The exclusive crushing container (made of the polypropylene) is an assembly type with double layer construction. The dust will not go out outside of the container during the crushing sieving in the device.
(outer measurement: H200mm, φ148mm)

Soil correction container with lid

The soil correction container can be used for sample keeping with lid.
(outer measurement: H73mm, φ141mm)

Product Information

Daiki Rika Kogyo Co., Ltd.

Head Office & Factory
212-8 Akagidai, Kounosu, Saitama,
365-0001 Japan
TEL:+81 48 568 2500 FAX:+81 48 568 2505
Branch Office
1-21 2-Chome, Nionohama, Otsu, Siga,
520-0801 Japan
TEL:+81 77 510 8550 FAX:+81 77 510 8555

<http://www.daiki.co.jp> e-mail:mbox@daiki.co.jp

DIK-2610
Dust shield Automatic mill and screen for soil RK-4 II